





Eclipta Prostrata Whole Plant – BP 2016

These chromatograms are provided for information only as an aid to analysts and intended as guidance for the interpretation and application of BP monographs.

Typical chromatogram for solution (1 & 2) in the Identification test for Eclipta Prostrata Whole Plant as published in BP 2016.




Typical chromatogram for solution (2) in the Assay test for Eclipta Prostrata Whole Plant as published in BP 2016.



The chromatogram shows two peaks, wedelolactone (60 mAU) at 16.0 minutes and coumestrol (50 mAU) at 27.5 minutes.



Typical chromatogram for solution (1) in the Assay test for *Eclipta Prostrata* Whole Plant as published in BP 2016.



The chromatogram shows one known peak, wedelolactone (43 mAU) at 15.5 minutes.