




Irbesartan Tablets – BP 2016


These chromatograms are provided for information only as an aid to analysts and intended as guidance for the interpretation and application of BP monographs.

Typical chromatogram for solution (5) in the Related Substances test for Irbesartan Tablets as published in BP 2016.



The chromatogram shows a large irbesartan peak (off scale at 50 mAU) at 5.8 minutes and one small peak due to impurity A (14 mAU) at 4.6 minutes.

Typical chromatogram for solution (1) in the Related Substances test for Irbesartan Tablets as published in BP 2016.



The chromatogram shows a large irbesartan peak (off scale at 25 mAU) at 5.8 minutes and one small peak due to impurity A (1 mAU) at 4.6 minutes.