




Testosterone Enantate Injection – BP 2016

This chromatogram is provided for information only as an aid to analysts and intended as guidance for the interpretation and application of BP monographs.

Typical chromatogram for solution (4) in the Related substances test for Testosterone Enantate Injection as published in BP 2016.



The chromatogram shows a large testosterone enantate peak (off scale at 40 mAU) at 21.5 minutes and two small peaks due to testosterone (37 mAU) at 1.5 minutes and testosterone caproate (2 mAU) at 15.0 minutes.